


Ifng

Mouse Interferon gamma ELISA Kit

Catalog No.	CDK345A CDK345B	Quantity:	1 x 96 tests 2 x 96 tests
Alternate Names:	IFN- γ , IFN-g, IFN-gamma		
Description:	The Cell Sciences™ IFN-g ELISA kit is a solid phase sandwich ELISA for the in-vitro qualitative and quantitative determination of murine IFN-g in supernatants, buffered solutions or serum and plasma samples. This assay will recognize both native and recombinant murine IFN-g.		
UniProt ID:	P01580		
Range:	31.25 - 1000 pg/ml		
Sensitivity:	15 pg/ml		
Sample Size:	100 μ l		
Incubation Time:	3 hours 45 minutes		
Precision:	Intra-assay CV 4.7% Inter-assay CV 5.2%		
Contents:	ELISA Kit includes pre-coated strip plate, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents and a detailed procedure.		
Storage & Stability:	Store entire kit at 2-8°C.		
Cross-Reactivity:	No cross reactivity with other murine cytokines.		



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

