

TNFRSF10D

Mouse Anti-Human TRAIL-R4 / CD264 (Clone B-P30) mAb, Azide Free

Catalog No.	CDM082A CDM082B	Quantity:	200 µg 500 µg
Alternate Names:	Tumor Necrosis Factor Receptor Superfamily Member 10D, Decoy Receptor 2, TNF-related apoptosis-inducing ligand receptor 4, TRAIL receptor 4, TRAIL receptor with a truncated death domain		
Description:	Tumor necrosis factor receptor superfamily member 10d is a member of the TNF-receptor superfamily. The receptor contains an extracellular TRAIL-binding domain, a transmembrane domain, and a truncated cytoplasmic death domain. The receptor does not induce apoptosis and has been shown to play an inhibitory role in TRAIL-induced cell apoptosis.		
UniProt ID:	Q9UBN6		
Gene ID:	8793		
Concentration:	1.0 mg/ml		
Specificity:	Recognizes TRAIL-R4, a 44.2 kDa protein		
Hybridoma:	Myeloma X63/AG.8653 x BALB/c spleen cells		
Immunogen:	Recombinant human TRAIL-R4/Fc chimera		
Isotype:	Mouse IgG1 κ		
Clone:	B-P30		
Formulation:	Sterile filtered PBS, carrier and preservative free.		
Applications:	Western Blot Immunoprecipitation		
Storage & Stability:	Store at 2-8°C for up to 1 year or in working aliquots at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com