


CD34

Mouse Anti-Human CD34 Class III (Clone B-F23) mAb

Catalog No.	CDM179	Quantity:	200 tests
Alternate Names:	CD34 antigen, Hematopoietic progenitor cell antigen CD34		
UniProt ID:	P28906		
Gene ID:	947		
Hybridoma:	Myeloma X63/AG.8653 x BALB/c spleen cells		
Specificity:	Recognizes the Gp105-120 antigen, a 105-120 kDa protein		
Host:	Mouse		
Immunogen:	KG-1a cell line		
Isotype:	IgG1		
Clone:	B-F23		
Concentration:	100 tests / 1 ml		
Formulation:	Liquid in PBS, 5% BSA, 0.1% sodium azide. Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Purification:	Ion exchange chromatography		
Applications:	Flow Cytometry Suitable for Immunohistochemistry		
Application Notes:	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood. The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Stable for 1 year at 2-8°C. Do Not Freeze.		



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com