


CD40LG

Mouse Anti-Human CD154 / CD40L (Clone B-B29) FITC mAb

Catalog No.	CDM184	Quantity:	100 tests
Alternate Names:	CD40 ligand, Tumor necrosis factor ligand superfamily member 5, TNFSF5, TNF-related activation protein, TRAP, T-cell gp39		
UniProt ID:	P29965		
Gene ID:	959		
Concentration:	100 tests/ml		
Specificity:	Recognises the CD40L, TRAP, T-BAM antigen, a 28/30/33 kDa protein		
Hybridoma:	Myeloma X63/AG.8653 x BALB/c spleen cells		
Immunogen:	CD40L transfected 3T6 cells		
Isotype:	Mouse IgG2a		
Clone:	B-B29		
Conjugate:	FITC		
Formulation:	PBS, 1% BSA, 0.1% sodium azide		
Purification:	Ion exchange chromatography		
Applications:	Suitable for flow cytometry		
Application Notes:	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood.		
Storage & Stability:	Stable at 2-8°C for 12 months. For longer storage, freeze aliquots. Avoid repeated freeze-thaw cycles.		

A typical staining pattern with the B-B29 monoclonal antibody of CD40L transfected 3T6 cell line.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com