

Mouse IgG1 Isotype Control (Clone B-Z1) mAb, Azide Free

Catalog No.	CDM377A CDM377B	Quantity:	200 µg 500 µg
Description:	The mouse monoclonal antibody is for use as an isotype control.		
Hybridoma:	Myeloma X63/AG.8653 x Balb/c spleen cells		
Concentration:	1 mg/ml		
Specificity:	Not specific for any known human antigens.		
Host:	Mouse		
Isotype:	IgG1κ		
Immunogen:	Rat monocytes and macrophages.		
Clone:	B-Z1		
Formulation:	Sterile-filtered PBS, treated to remove endotoxins. Carrier and preservative free.		
Applications:	Flow Cytometry, IHC, Negative control for Functional Assays.		
Storage & Stability:	Stable at 2-8°C for 1 year. For long term storage, freeze in aliquots at -20 to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

