

Mouse IgG2a Isotype Control (Clone B-Z2) mAb, Azide Free

Catalog No.	CDM383A CDM383B	Quantity:	200 µg 500 µg
Description:	The mouse monoclonal antibody is for use as an isotype control.		
Concentration:	1 mg/ml		
Specificity:	Negative control for biological test.		
Hybridoma:	Myeloma X63/AG.8653 x Balb/c spleen cells		
Isotype:	Mouse IgG2a κ		
Immunogen:	Rat monocytes and macrophages.		
Clone:	B-Z2		
Formulation:	Sterile-filtered PBS, carrier and preservative free.		
Cross-Reactivity:	Not specific for any known human antigens		
Applications:	Flow Cytometry		
Storage & Stability:	Store at 2-8°C for up to 1 year or in working aliquots at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

