


IL1RAP

Mouse Anti-Human IL-1RAcP / IL-1R3 (Clone B-L43) mAb, Azide Free

Catalog No.	CDM524A CDM524B	Quantity:	200 µg 500 µg
Alternate Names:	Interleukin-1 receptor accessory protein, IL-1RAcP, Interleukin-1 receptor 3, IL-1R-3, IL-1R3		
Description:	<p>The pro-inflammatory cytokine IL-1 induced cellular response requires two subunits of its receptor, IL-1 receptor I(IL-1RI) and IL-1 receptor accessory protein (IL-1RAcP). The IL-1 Receptor Accessory Protein (IL-1RAcP) forms a complex with IL-1RI in response to IL-1 treatment. The IL-1 receptor-associated kinase (IRAK), which mediates activation of NF-κB inducing kinase (NIK) and of NF-κB, recruits to the IL-1R complex through IL-1 RAcP. IL-1 activation of stress-activated protein kinase and of acid sphingomyelinase also requires IL-1RAcP. Like IL-1RI, IL-1RAcP subunit is essential for IL-1 mediated cellular response. IL-1RAcP is expressed in many tissues.</p>		
UniProt ID:	Q9NPH3		
Gene ID:	3556		
Hybridoma:	Myeloma X63/AG.8653 x Balb/c spleen cells		
Specificity:	Recognizes native and recombinant human IL1RAP		
Isotype:	Mouse IgG1 kappa		
Immunogen:	Recombinant human IL1RAcP		
Clone:	B-L43		
Concentration:	1.0 mg/ml		
Formulation:	Sterile-filtered PBS, carrier and preservative free.		
Applications:	ELISA Flow cytometry		
Storage & Stability:	Stable at 2-8°C for 12 months. For longer storage, freeze aliquots at -20 to -80°C. Avoid repeated freeze-thaw cycles.		



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com