

Human IL-17A Matched Antibody Pair for ELISPOT

Catalog No: CK1001 **Size:** 10 plate format

Coating antibodies (2 vials)

Product: Monoclonal antibody to human interleukin 17A (IL-17A)
Isotype: Mouse IgG2b
Production: *In vitro* using serum free medium
Purification: Ion exchange chromatography
Contents: Each vial contains sufficient material for coating of five 96-well ELISPOT plates
Buffer: Prior to lyophilization: PBS + 125 mM trehalose
Application: Coating antibody in an ELISPOT system

Detector antibodies (2 vials)

Product: Biotinylated polyclonal antibody to human interleukin 17A (IL-17A)
Isotype: Mouse IgG2b
Production: *In vitro* using serum free medium
Purification: Ion exchange chromatography
Labeling: With biotin (biotinamidohexanoic acid N-hydroxysuccinimide ester)
Contents: Each vial contains sufficient material for five 96-well ELISPOT plates
Buffer: Prior to lyophilization: PBS + 1% BSA + 125 mM trehalose
Application: Detection antibody in an ELISPOT system

General *Please read our ELISPOT antibody pair Technical Data Sheet for guidelines and procedures.*

Specificity: Validated for detecting natural human IL-17A
Sterility: Membrane filtered (0.2 µm)
Content: Lyophilized
Stability: The lyophilized products are stable for more than one year at 2-8°C. After reconstitution, the antibodies are stable for several months at 2-8°C (if kept sterile) or for minimal one year at -20°C.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences, Inc.
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888 769-1246
Phone: 978 572-1070
Fax: 978 992-0298

E-mail: info@cellsciences.com
Web Site: www.cellsciences.com