

## Monkey G-CSF Matched Antibody Pair for ELISPOT

**Catalog No:** CK1003                      **Size:** 10-plate format

### Coating antibodies (2 vials)

**Product:** Monoclonal antibody to monkey granulocyte-colony stimulating factor (G-CSF)  
**Isotype:** Mouse IgG1  
**Production:** *In vitro* using serum free medium  
**Purification:** Ion exchange chromatography  
**Contents:** Each vial contains sufficient material for coating of five 96-well ELISPOT plates  
**Buffer:** Prior to lyophilization: PBS + 125 mM trehalose  
**Application:** Coating antibody in an ELISPOT system

### Detector antibodies (2 vials)

**Product:** Biotinylated polyclonal antibody to monkey granulocyte-colony stimulating factor (G-CSF)  
**Isotype:** Rabbit IgG  
**Purification:** Ammonium sulphate precipitation and protein A-affinity chromatography  
**Labeling:** With biotin (biotinamido-hexanoic acid N-hydroxysuccinimide ester)  
**Contents:** Each vial contains sufficient material for five 96-well ELISPOT plates  
**Buffer:** Prior to lyophilization: PBS + 1% BSA + 125 mM trehalose  
**Application:** Detection antibody in an ELISPOT system

**General**                      *Please read our ELISPOT antibody pair Technical Data Sheet for guidelines and procedures.*

**Specificity:** Validated for detecting monkey G-CSF  
**Sterility:** Membrane filtered (0.2 µm)  
**Content:** Lyophilized  
**Stability:** The lyophilized products are stable for more than one year at 2-8°C. After reconstitution, the antibodies are stable for several months at 2-8°C (if kept sterile) or for one year at -20°C.

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**



**Cell Sciences, Inc.**  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888 769-1246  
Phone: 978 572-1070  
Fax: 978 992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Web Site: [www.cellsciences.com](http://www.cellsciences.com)