

Monkey IFN-gamma Matched Antibody Pair for ELISPOT

Catalog No: CK1004 **Size:** 10-plate format

Coating antibodies (2 vials)

Product: Monoclonal antibody to monkey interferon gamma (IFN- γ)
Isotype: Mouse IgG1
Production: *In vitro* using serum free medium
Purification: Ion exchange chromatography
Contents: Each vial contains sufficient material for coating of five 96-well ELISPOT plates
Buffer: Prior to lyophilization: PBS + 125 mM trehalose
Application: Coating antibody in an ELISPOT system

Detector antibodies (2 vials)

Product: Biotinylated Monoclonal antibody to monkey interferon gamma (IFN- γ)
Isotype: Mouse IgG1
Production: *In vitro* using serum free medium
Purification: Ion exchange chromatography
Labeling: With biotin (biotinamido hexanoic acid N-hydroxysuccinimide ester)
Contents: Each vial contains sufficient material for five 96-well ELISPOT plates
Buffer: Prior to lyophilization: PBS + 1% BSA + 125 mM trehalose
Application: Detection antibody in an ELISPOT system

General *Please read our ELISPOT antibody pair Technical Data Sheet for guidelines and procedures.*

Specificity: Validated for detecting rhesus macaque, cynomolgus monkey, pig-tailed macaque, baboon, African green monkey IFN- γ
Sterility: Membrane filtered (0.2 μ m)
Content: Lyophilized
Stability: The lyophilized products are stable for more than one year at 2-8°C. After reconstitution, the antibodies are stable for several months at 2-8°C (if kept sterile) or for one year at -20°C.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences, Inc.
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888 769-1246
Phone: 978 572-1070
Fax: 978 992-0298

E-mail: info@cellsciences.com
Web Site: www.cellsciences.com