

## IFNGR2

### Mouse Anti-Human IFN-gamma Receptor 2 (Clone MMHGR-2) mAb

<b>Catalog No.</b>	CM2027	<b>Quantity:</b>	250 µg
<b>Alternate Names:</b>	AF-1, IFGR2, IFNGT1, IMD28		
<b>Description:</b>	Mouse Anti-Human IFN-gamma Receptor 2 Clone MMHGR-2 mAb. This is a non-neutralizing antibody.		
<b>Concentration:</b>	> 1 mg/ml, lot specific		
<b>Gene ID:</b>	3460		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Specificity:</b>	Human IFN-gamma Receptor Chain 2. Binds to but does not neutralize human interferon gamma receptor.		
<b>Host:</b>	Mouse		
<b>Immunogen:</b>	Human IFN-gamma Receptor Chain 2		
<b>Isotype:</b>	IgG1		
<b>Clone:</b>	MMHGR-2		
<b>Endotoxin:</b>	< 1 EU/µg		
<b>Formulation:</b>	Liquid in PBS containing 0.1% BSA.		
<b>Applications:</b>	Western Blot: working dilution of 1:2,000 Immunoprecipitation: working dilution of 1:100 Immunohistochemistry: working dilution of 1:100 Flow Cytometry: working dilution of 1:50 The optimal concentration should be determined by the user for each specific application.		
<b>Storage &amp; Stability:</b>	<b>Centrifuge vial prior to opening.</b> This product should be kept at -20°C or below for retention of full activity. When thawing, the contents of the tube should be apportioned in separate tubes. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences**<sup>®</sup>  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)