

## Mouse Anti-HA Clone 4B2 mAb

<b>Catalog No.</b>	CMAB1004A	<b>Quantity:</b>	100 µg
	CMAB1004B		200 µg
	CMAB1004C		1 mg

**Description:** Epitope tags are useful for the labeling and detection of proteins using immunoblotting, immunoprecipitation, and immunostaining techniques. The HA tag is derived from an epitope of the *Influenza* Hemagglutinin (HA) protein, which has been extensively used as a general epitope tag in expression vectors.

**Concentration:** 0.5 mg/mL

**Host:** Mouse

**Isotype:** IgG1

**Clone:** 4B2

**Formulation:** Liquid in PBS + 50% Glycerol

**Applications:** WB, IP

**Application Notes:** Recommended dilution for WB: 1.0 µg/mL  
The optimal concentration should be determined by the user for each specific application.

**Storage & Stability:** Store at -20°C. **Avoid repeated freeze-thaw cycles.**

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

