



Mouse Anti-Rabbit IgG (Clone 2A9) HRP mAb

Catalog No.	CMB006	Quantity:	1.0 mL
Description:	Monoclonal mouse anti-Rabbit IgG (Clone 2A9) HRP-labeled is validated for ELISA and Western blot.		
Specificity:	Rabbit IgG		
Immunogen:	Recombinant rabbit IgG Fc fragment		
Isotype:	Mouse (BALB/c) IgG ₁ κ		
Clone:	2A9		
Conjugate:	HRP (Horseradish Peroxidase)		
Formulation:	50% glycerol/50% PBS, pH 7.4. No preservative added.		
Isotype Control	Mouse IgG ₁ -HRP		
Cross-Reactivity:	Minimal reactivity to human, mouse, rat, goat, bovine, horse, donkey, sheep, guinea pig, and chicken serum.		
Applications:	ELISA, Immunohistochemistry, Western blot		
Application Notes:	Recommended dilution for ELISA is 1:4,000 - 1:8,000		
Storage & Stability:	Store at 2-8°C or long-term at -20°C to -80°C .		

ELISA plate was coated with purified rabbit IgG and IgM. Immunoglobulins were detected with serially diluted Mouse Anti-Rabbit IgG-HRP.

Rabbit immunoglobulins above were resolved by electrophoresis under reducing conditions, transferred to PVDF membrane, and probed with Mouse Anti-Rabbit IgG-HRP. Proteins were visualized using chemiluminescent detection.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.