


Mouse Anti-Human IgM (Clone UHB) mAb, HRP

Catalog No.	CMH000	Quantity:	1.0 mL
Description:	Monoclonal antibody Mouse Anti-Human IgM Clone UHB, Horseradish peroxidase-labeled		
Specificity:	Human IgM (900 kDa)		
Host:	Mouse		
Isotype:	Mouse IgG ₃ k		
Clone:	UHB		
Conjugate:	HRP (Horseradish Peroxidase)		
Formulation:	Liquid stock solution in 50% glycerol/50% PBS, pH 7.4. No preservative added.		
Isotype Control	Mouse IgG ₃		
Applications:	ELISA, Western Blot		
Application Notes:	ELISA: suggested dilution range 1:4,000 – 1:8,000. The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Store unopened at 2-8°C for up to 1 year. Upon opening, it is recommended to store in working aliquots at -20°C to -80°C. Avoid repeated freeze/thaw cycles.		

ELISA plate was coated with purified human IgM, IgG, and IgA. Immunoglobulins were detected with serially diluted Mouse Anti-Human IgM-HRP.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

