

Mouse Anti-Human IgG₂ Fd Clone HP6014 Biotin mAb

Catalog No.	CMH024	Quantity:	500 µg
Concentration:	500 µg/mL		
Specificity:	Human IgG ₂ Fd.		
Immunogen:	Human IgG ₂ myeloma protein.		
Isotype:	Mouse (BALB/c) IgG ₁ k		
Clone:	HP6014		
Conjugate:	Biotin		
Formulation:	Liquid volume of 1.0 mL in a buffer of PBS/NaN ₃ .		
Isotype Control	Mouse IgG ₁ -Biotin		
Applications:	Applications for relevant formats of this clone include ELISA, FLISA, Flow Cytometry, Immunohistochemistry, Western Blot, and Immunoprecipitation.		
Application Notes:	Working dilutions for ELISA 1:1,000 – 1:2,000 , for FLISA with FTIC or AF488 conjugates 1:200 - 1:400. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C. For long term storage at -20°C, aliquot into smaller volumes to avoid multiple freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

