


Mouse Anti-Human Lambda Clone JDC-12 HRP mAb

Catalog No.	CMH029	Quantity:	1.0 mL
Description:	Mouse Anti-Human Lambda Clone JDC-12 HRP mAb		
Specificity:	Human/Rhesus Lambda		
Isotype:	Mouse IgG1κ		
Clone:	JDC-12		
Conjugate:	HRP (Horseradish Peroxidase)		
Formulation:	Liquid volume of 1.0 mL in a stock solution of 50% glycerol/50% PBS, pH 7.4. No preservative added.		
Isotype Control	Mouse IgG1-HRP		
Applications:	Applications for relevant formats of this clone include ELISA and Western Blot.		
Application Notes:	Working dilutions for ELISA 1:1,000 – 1:2,000. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C or long-term at -20°C.		

ELISA plate was coated with purified human IgGκ, IgMκ, IgAκ, IgGλ, IgMλ, and IgAλ. Immunoglobulins were detected with serially diluted Mouse Anti-Human Lambda-HRP.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

