


Mouse Anti-Human IgG₃ Hinge (Clone HP6050) mAb, HRP

Catalog No.	CMH033	Quantity:	1.0 mL
Description:	Mouse Anti-Human IgG ₃ Hinge (Clone HP6050) HRP mAb may be used for ELISA, Western blot and Immunohistochemistry applications.		
Specificity:	Human IgG ₃ Hinge, Mr 170 kDa		
Immunogen:	Human IgG ₃ myeloma protein		
Isotype:	Mouse (BALB/c) IgG ₁ k		
Clone:	HP6050		
Conjugate:	HRP (Horseradish Peroxidase)		
Formulation:	50% glycerol/50% PBS, pH 7.4. No preservative added.		
Isotype Control	Mouse IgG ₁ -HRP		
Application Notes:	ELISA: suggested dilution 1:4,000 – 1:10,000 The optimum working dilution should be determined by the user for their specific application.		
Storage & Stability:	Store at 2-8°C for immediate use or in working aliquots at -20°C to -80°C for up to 1 year.		

ELISA plate was coated with purified human IgG1, IgG2, IgG3, and IgG4. Immunoglobulins were detected with serially diluted Mouse Anti-Human IgG3 Hinge-HRP.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

