


Mouse Anti-Human IgG₃ Hinge (Clone HP6050) mAb, Biotin

Catalog No.	CMH034	Quantity:	500 µg
Description:	Monoclonal mouse Anti-Human IgG ₃ Hinge region (Clone HP6050) with Biotin-label		
Concentration:	0.5 mg/ml		
Volume:	1.0 ml		
Specificity:	Human IgG ₃ Hinge (Mr 170 kDa)		
Immunogen:	Human IgG3 myeloma protein		
Isotype:	Mouse (BALB/c) IgG ₁ κ		
Clone:	HP6050		
Conjugate:	Biotin		
Formulation:	Liquid PBS/NaN ₃		
Isotype Control	Mouse IgG ₁ -Biotin		
Applications:	ELISA, Western blot		
Application Notes:	ELISA 1:5,000 – 1:10,000 working dilution. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C.		

ELISA plate was coated with purified human IgG1, IgG2, IgG3, and IgG4. Immunoglobulins were detected with serially diluted Mouse Anti-Human IgG3 Hinge-BIOT followed by Streptavidin-HRP.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

