

## Mouse IgG1 Isotype Control (Clone 15H6) mAb

<b>Catalog No.</b>	CMH050	<b>Quantity:</b>	500 µg
<b>Description:</b>	Low Endotoxin/Azide Free mouse monoclonal antibody useful as an mouse IgG1 isotype control.		
<b>Concentration:</b>	0.5 mg/ml		
<b>Conjugate:</b>	Unconjugated		
<b>Specificity:</b>	Not specific for any known human antigens.		
<b>Host:</b>	Mouse		
<b>Isotype:</b>	IgG1		
<b>Clone:</b>	15H6		
<b>Formulation:</b>	Sterile filtered PBS, pH 7.4 without preservative		
<b>Endotoxin:</b>	< 0.1 EU/µg		
<b>Storage &amp; Stability:</b>	Stable at 2-8°C for 1 year. Handle under aseptic conditions		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

