

IFNG

Mouse Anti-human IFN-gamma (Clone MD-2) mAb

Catalog No.	CMI018	Quantity:	0.5 mg
Alternate Names:	IFN-gamma, Immune interferon		
Description:	Mouse Anti-Human IFN-gamma Clone MD-2, monoclonal antibody produced <i>in vitro</i> using serum free medium.		
UniProt ID:	P01579		
Gene ID:	3458		
Specificity:	Binds with high affinity to both natural and recombinant human IFN-gamma but does not neutralize the antiviral bioactivity of IFN-gamma <i>in vitro</i> . The antibody exerts no reactivity towards mouse and rat IFN-gamma and human IFN-alpha and -beta.		
Source:	Produced <i>in vitro</i> using serum free medium.		
Host:	Mouse		
Isotype:	IgG ₁		
Clone:	MD-2		
Purification:	Ion exchange chromatography		
Formulation:	Lyophilized from sterile-filtered PBS containing 125 mM Trehalose		
Reconstitution:	Centrifuge vial briefly before opening. Reconstitute with 0.5 ml sterile distilled water, containing 0.02% sodium azide to prevent bacterial growth (recommended). Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Quantitation:	$E^{0.1\%}_{280\text{ nm}} = 1.4$ for a 1 mg/ml solution.		
Applications:	Coating antibody in ELISA Immunohistochemistry Western blot analysis		
Storage & Stability:	Lyophilized product is stable for at least one year at 2-8°C. After reconstitution, the contents can be safely stored at 2-8°C for one month or for one year in working aliquots at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com