

CCL2

Mouse Anti-Human MCP-1 / CCL2 (Clone B-Z21) mAb

Catalog No.	CMI036	Quantity:	0.5 mg
Alternate Names:	C-C motif chemokine 2, CCL2, Monocyte chemotactic protein 1, MCP-1, Monocyte chemotactic and activating factor, MCAF, Monocyte secretory protein JE, Small-inducible cytokine A2		
Description:	This monoclonal antibody (clone B-Z21) to human MCP-1 is useful for ELISA		
UniProt ID:	P13500		
Immunogen:	Recombinant human MCP-1		
Source:	Produced <i>in vitro</i> using serum free medium.		
Specificity:	Recognizes native and recombinant human MCP-1		
Isotype:	Mouse IgG1 kappa		
Clone:	B-Z21		
Quantitation:	$E^{0.1\%}_{280\text{ nm}} = 1.4$ for a 1 mg/ml solution.		
Formulation:	Lyophilized from sterile-filtered PBS containing 125 mM Trehalose		
Purification:	Ion exchange chromatography.		
Reconstitution:	Centrifuge vial briefly before opening. Reconstitute with 0.5 ml sterile distilled water, containing 0.02% sodium azide to prevent bacterial growth (recommended).		
Application:	ELISA		
Storage & Stability:	Lyophilized product is stable for at least one year at 2-8°C. After reconstitution, the contents can be safely stored at 2-8°C for one month or for one year in working aliquots at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com