

## CCL2

### Mouse Anti-Human MCP-1 / CCL2 (Clone B-K41) mAb

<b>Catalog No.</b>	CMI037	<b>Quantity:</b>	0.5 mg
<b>Alternate Names:</b>	C-C motif chemokine 2, CCL2, Monocyte chemotactic protein 1, MCP-1, Monocyte chemotactic and activating factor, MCAF, Monocyte secretory protein JE, Small-inducible cytokine A2		
<b>Description:</b>	This monoclonal antibody (clone B-K41) to human MCP-1 is useful for ELISA		
<b>UniProt ID:</b>	P13500		
<b>Immunogen:</b>	Recombinant human MCP-1		
<b>Source:</b>	Produced <i>in vitro</i> using serum free medium.		
<b>Specificity:</b>	Recognizes native and recombinant human MCP-1		
<b>Isotype:</b>	Mouse IgG1 kappa		
<b>Clone:</b>	B-K41		
<b>Quantitation:</b>	$E^{0.1\%}_{280\text{ nm}} = 1.4$ for a 1 mg/ml solution.		
<b>Formulation:</b>	Lyophilized from sterile-filtered PBS containing 125 mM Trehalose		
<b>Purification:</b>	Ion exchange chromatography.		
<b>Reconstitution:</b>	<b>Centrifuge vial briefly before opening.</b> Reconstitute with 0.5 ml sterile distilled water, containing 0.02% sodium azide to prevent bacterial growth (recommended).		
<b>Application:</b>	ELISA		
<b>Storage &amp; Stability:</b>	Lyophilized product is stable for at least one year at 2-8°C. After reconstitution, the contents can be safely stored at 2-8°C for one month or for one year in working aliquots at -20°C to -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences**<sup>®</sup>  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)