

Mouse Anti-Monkey IgG (Clone SB108a) HRP mAb

Catalog No.	CMK000	Quantity:	1.0 mL
Description:	Mouse Anti-Monkey IgG, Clone SB108a, horseradish peroxidase-labeled monoclonal antibody for ELISA		
Specificity:	Rhesus and Cynomolgus IgG		
Immunogen:	Rhesus IgG		
Isotype:	Mouse (BALB/c) IgG1k		
Clone:	SB108a		
Conjugate:	HRP (Horseradish Peroxidase)		
Formulation:	50% glycerol/50% PBS, pH 7.4 No preservative added.		
Isotype Control	Mouse IgG1-HRP		
Cross-Reactivity:	Minimal reactivity to human, mouse, rat, hamster, guinea pig, rabbit, goat, sheep, donkey, bovine, horse, porcine, feline, and chicken serum		
Applications:	ELISA: suggested starting dilution 1:4,000-1:8,000. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C or long-term at -20°C.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

