

Mouse Anti-Monkey IgG₁ Clone SB108a Biotin mAb

Catalog No.	CMK001	Quantity:	1.0 mL
Description:	Mouse Anti-Monkey IgG ₁ Clone SB108a Biotin mAb		
Concentration:	500 µg/mL		
Specificity:	Rhesus and cynomolgus IgG		
Immunogen:	Rhesus IgG		
Isotype:	Mouse (BALB/c) IgG1κ		
Clone:	SB108a		
Conjugate:	HRP (Horseradish Peroxidase)		
Formulation:	Liquid volume of 1.0mL in a buffer of PBS/NaN ₃		
Isotype Control	Mouse IgG ₁ -Biotin		
Cross-Reactivity:	Minimal reactivity to human, mouse, rat, hamster, guinea pig, rabbit, goat, sheep, donkey, bovine, horse, porcine, feline, and chicken serum		
Applications:	Applications for relevant formats of this clone include Flow Cytometry, ELISA and FLISA.		
Application Notes:	Working dilutions; ELISA BIOT conjugate 1:5,000 - 1:10,000. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C. For long term storage at -20°C, aliquot into smaller volumes to avoid multiple freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

