

## Mouse Anti-Monkey IgG (Clone SB108a) mAb

<b>Catalog No.</b>	CMK002	<b>Quantity:</b>	0.5 mg
<b>Description:</b>	Mouse Anti-Monkey IgG (Clone SB108a) mAb reacts with Rhesus macaque and Cynomolgus monkey IgG. Minimal reactivity with Human IgG.		
<b>Immunogen:</b>	Rhesus macaque IgG		
<b>Specificity:</b>	Rhesus macaque and cynomolgus monkey IgG		
<b>Cross-Reactivity:</b>	Minimal reactivity to human, mouse, rat, hamster, guinea pig, rabbit, goat, sheep, donkey, bovine, horse, porcine, feline, and chicken serum		
<b>Concentration:</b>	0.5 mg/ml		
<b>Volume:</b>	1.0 ml		
<b>Isotype:</b>	Mouse (BALB/c) IgG1 $\kappa$		
<b>Clone:</b>	SB108a		
<b>Formulation:</b>	Borate buffered saline, pH 8.2, without carrier or preservative.		
<b>Applications:</b>	ELISA, Flow cytometry, Western blot		
<b>Storage &amp; Stability:</b>	Store at 2-8°C. For long term storage at -20°C, aliquot into smaller volumes to avoid multiple freeze/thaw cycles.		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

