


Rat Anti-Mouse IgG2a Biotinylated mAb

Catalog No.	CMM026	Quantity:	500 µg
Description:	Rat Anti-Mouse IgG _{2a} Clone SB84a Biotin mAb		
Concentration:	500 µg/mL		
Specificity:	Mouse IgG _{2a}		
Immunogen:	Mouse IgG _{2a} hybridoma		
Isotype:	Rat (Sprague Dawley) IgG _{1κ}		
Clone:	SB84a		
Conjugate:	Biotin		
Formulation:	PBS containing NaN ₃		
Isotype Control	Rat IgG ₁ -BIOT		

Applications: **Flow Cytometry:** suggested concentration < 1 µg/10⁶ cells
ELISA: suggested dilution 1:5,000 - 1:20,000
 Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.

Storage & Stability: Store at 2-8°C. For long term storage at -20°C, aliquot into smaller volumes to avoid multiple freeze/thaw cycles.

ELISA plate was coated with purified mouse IgG1, IgG2a, IgG2b, IgG3, IgM, and IgA. Immunoglobulins were detected with serially diluted Rat Anti-Mouse IgG2a-Biotin followed by Streptavidin-HRP (Cat. No. CPO016).



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.