


Mouse Anti-Rat IgG2b Clone 2B10A8 PE mAb

Catalog No.	CMR021	Quantity:	100 µg
Description:	Mouse Anti-Rat IgG2b Clone 2B10A8 PE mAb		
Concentration:	100 µg/mL		
Specificity:	Rat IgG _{2b}		
Immunogen:	Rat IgG _{2b} hybridoma		
Isotype:	Mouse (BALB/c) IgG _{2b} K		
Clone:	2B10A8		
Conjugate:	PE (R-phycoerythrin)		
Formulation:	Liquid volume of 1.0 mL in PBS/NaN ₃ and a stabilizing agent		
Isotype Control	Mouse IgG ₁ -PE		
Applications:	Applications for relevant formats of this clone include ELISA, FLISA, Flow Cytometry, Immunohistochemistry, Western Blot.		
Application Notes:	Working dilutions; FLISA ≤ 1 µg/ml, for Flow Cytometry ≤ 0.1 µg/10 ⁶ cells. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C. DO NOT FREEZE		

FLISA plate was coated with purified rat IgG1, IgG2a, IgG2b, IgG2c, and IgM. Immunoglobulins were detected with serially diluted Mouse Anti-Rat IgG2b-PE.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

