

Mouse Anti-human VEGFD mAb

Catalog No.	CMV116	Quantity:	100 µg
Description:	Mouse Anti-human VEGFD monoclonal antibody		
Alternate Names:	FIGF, Vascular Endothelial Growth Factor D, VEGF-D		
Specificity:	Recognizes human VEGFD		
Host:	Mouse		
Immunogen:	Human VEGFD recombinant protein		
Isotype:	IgG2		
Formulation:	Lyophilized from a 0.2 µm sterile filtered solution in PBS		
Purification:	Protein G affinity chromatography		
Reconstitution:	Centrifuge vial prior to opening. Reconstitute with sterile distilled water.		
Cross-Reactivity:	No cross-reactivity is shown with VEGFA and VEGFB.		
Applications	Western Blot Immunohistochemistry		
	The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Lyophilized antibody is stable at -80°C. Reconstituted antibody can also be aliquotted and stored frozen at -80°C for at least six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.