

VTN

Mouse Anti-Human Vitronectin Clone 4A132 mAb

Catalog No.	CMV135B CMV135C	Quantity:	1 mg 10 mg
Alternate Names:	V75, VN, VNT		
Gene ID:	7448		
Description:	Mouse Anti-Human Vitronectin Clone 4A132 monoclonal antibody		
Concentration:	10 mg/ml		
Specificity:	Human vitronectin		
Host:	Mouse		
Isotype:	IgG		
Clone:	4A132		
Formulation:	Liquid in 0.05 M sodium phosphate + 0.1 M NaCl + 1 mM EDTa, pH 6.6		
Purification:	Protein A purification		
Cross-Reactivity:	Does not cross-react with mouse.		
Applications:	Western Blot		
Application Notes:	Recognizes reduced and non-reduced monomeric human vitronectin in Western Blot. The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Store in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
480 Neponset Street
Bldg 12A
Canton, MA 02021

Toll Free: 888-769-1246
Phone: 781-828-0610
Fax: 781-828-0542

E-mail: info@cellsciences.com
Website: www.cellsciences.com