

Normal Cynomolgus Monkey Serum, Pooled, Off the Clot

Catalog No.	CNK001A	Quantity:	50 ml
	CNK001B		100 ml

Description: Normal Cynomolgus monkey serum is made from whole blood collected from normal, healthy cynomolgus monkeys born in the USA. The serum is removed from spontaneously-clotted whole blood which is centrifuged to separate the serum from the clot. Sera from at least 3 animals are pooled for this product.
Please inquire about single animal serum.

Volume: 50 mL

Appearance Yellow to orange liquid

Source: *Macaca fascicularis* whole blood, from at least 3 animals.

Storage & Stability: Store in working aliquots at -80 °C. Avoid repeated freeze-thaw cycles.

Origin: USA

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

