

PLG

Native Cynomolgus Monkey Plasminogen

Catalog No.	CNK002A CNK002B	Quantity:	1 mg 10 mg
Description:	Native Cynomolgus Macaque Monkey (<i>Macaca fascicularis</i>) Plasminogen prepared from frozen plasma by immobilized lysine chromatography.		
Concentration:	10.0 mg/ml		
Gene ID:	102145643		
Source:	Cyno monkey plasma		
Molecular Weight:	92 kDa		
Formulation:	Frozen liquid in 0.1 M HEPES + 0.1 M NaCl, pH 7.4.		
Purity:	>95% by SDS-PAGE analysis		
Storage & Stability:	Store at -80°C. Stable for 3 years from delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

