

Normal Cynomolgus Monkey Serum, Single Donor, Off the Clot

Catalog No.	CNK003A	Quantity:	1.0 ml
	CNK003B		10 ml

Description: Normal Cynomolgus monkey serum is made from whole blood collected from normal, healthy, adult cynomolgus monkeys born in the USA. The serum is removed from spontaneously-clotted whole blood which is centrifuged to separate the serum from the clot. Serum from a single, unique animal is collected for this product.

Appearance Yellow to orange liquid

Source: *Macaca fascicularis* whole blood

Storage & Stability: Store in working aliquots at -80 °C. Avoid repeated freeze-thaw cycles.

Origin: USA

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.