

Recombinant Streptavidin

Catalog No.	CPO013	Quantity:	0.5 mg
Description:	Purified streptavidin		
Concentration:	0.5 mg/ml		
UniProt ID:	P22629		
Source:	Recombinant <i>E. coli</i> -derived		
Specificity:	Biotin		
Conjugate:	Unconjugated		
Formulation:	Borate-buffered saline, pH 8.2. No preservatives or amine-containing buffer salts added.		
Applications:	ELISA		
Application Notes:	For coating ELISA plates, recommended dilution approximately 5 µg/mL. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

