

Streptavidin-Alkaline Phosphatase

Catalog No.	CPO015	Quantity:	1.0 mL
Description:	Streptavidin-AP		
Source:	Recombinant <i>E. coli</i> -derived		
Specificity:	Biotin		
Conjugate:	Alkaline Phosphatase, AP		
Formulation:	50 mM Tris, 1 mM MgCl ₂ , 50% glycerol, 8.0. Contains NaN ₃ as preservative.		
Applications:	ELISA		
Application Notes:	Working dilution; ELISA 1:2,000 - 1:4,000. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at 2-8°C or long term at -20 °C.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

