

Protein A

Catalog No.	CPO019	Quantity:	1.0 mg
Description:	Protein A		
Specificity:	Protein A binds the Fc portion of human IgG ₁ , IgG ₂ , IgG ₃ (weakly), IgG ₄ and mouse IgG ₁ (weakly), IgG _{2a} , IgG _{2b} and IgG ₃ ; Protein A also binds monkey, rabbit, porcine, guinea pig, canine and feline IgG.		
Formulation:	Lyophilized. No BSA, preservatives or amine-containing buffer salts have been added.		
Applications:	ELISA, Affinity Chromatography, Immunoprecipitation.		
Application Notes:	Working dilutions for ELISA 1:4,000 - 1:8,000. Since applications vary, you should determine the optimum working dilution for the product that is appropriate for your specific need.		
Storage & Stability:	Store at -20 °C for one year or reconstituted in buffered saline for one month at 2-8 °C .		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

