

TNNC1, TNNI3, TNNT2

Native Human Cardiac Troponin C-I-T Complex

Catalog No.	CRA146A CRA146B CRA146C CRA146D	Quantity:	25 µg 50 µg 1.0 mg 100 µg
Alternate Names:	Troponin C, slow skeletal and cardiac muscles, TN-C, Troponin I, cardiac muscle, Cardiac troponin I, Troponin T, cardiac muscle, TnTc, Cardiac muscle troponin T, cTnT		
Description:	Troponin Complex is a heteromeric protein playing an important role in the regulation of skeletal and cardiac muscle contraction. It consists of three subunits, Troponin I, Troponin T and Troponin C. Each subunit is responsible for part of Troponin Complex function. E.g. Troponin I inhibits ATP-ase activity of acto-myosin. Troponin T and Troponin I are presented in cardiac muscles in different forms than in skeletal muscles.		
Uniprot ID:	P63316, P19429, P45379		
Source:	Human cardiac tissue		
Molecular Weight:	approximately 75 kDa		
Formulation:	Sterile-filtered 25 mM Tris buffer, 0.15 M NaCl, 0.09% sodium azide pH 7.5		
Purity:	> 90% as determined by SDS-PAGE analysis (Coomassie)		
Applications:	Western Blot Immunoassays		
Storage & Stability:	Store at 2-8°C up to 2 weeks. For longer term, store aliquots at -20°C to -80°C, containing 0.1% HSA or BSA. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.