

Recombinant Mouse Alpha 2 Antiplasmin with Flag Tag

Catalog No.	CRA159B CRA159C	Quantity:	100 µg 1 mg
Description:	Recombinant Mouse Alpha 2 Antiplasmin with FLAG tag		
Concentration:	5.0 mg/ml		
Source:	Insect cells		
Molecular Weight:	65.0 kDa		
Formulation:	Liquid in 0.1 M Tris-HCl + 0.15 M NaCl, pH 7.4		
Purification:	FLAG affinity		
Storage & Stability:	Store in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

