

Recombinant Human AKT3/Protein Kinase B Gamma GST Active

Catalog No.	CRA215A CRA215B	Quantity:	5 µg 10 µg
Alternate Names:	Protein kinase B-gamma, PKBG, PKB-gamma, PRKBG, RAC-gamma, STK-2		
Description:	Recombinant full-length human AKT3 containing N-terminal GST tag.		
Concentration:	0.1 mg/ml		
GeneID:	10000		
Source:	Sf9 insect cells		
Molecular Weight:	~85 kDa		
Formulation:	Liquid in 50 mM Tris-HCl, pH 7.5 + 150 mM NaCl + 0.25 mM DTT + 0.1 mM EGTA + 0.1 mM EDTA + 0.1 mM PMSF + 25% glycerol		
Purity:	> 90%		
Specific Activity:	206 nmol/min/mg		
Storage & Stability:	Stable for 1 year in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

