

Recombinant Human CHEK1 GST Active

Catalog No.	CRC109A CRC109B	Quantity:	5 µg 10 µg
Alternate Names:	Checkpoint Kinase 1 homolog, CHK1		
Description:	Recombinant full length human CHEK1 containing N-terminal GST tag.		
Concentration:	0.1 mg/ml		
GenelD:	1111		
Protein Accession No:	NM_001274		
Source:	Sf9 insect cells		
Molecular Weight:	~82 kDa		
Formulation:	Liquid in 50 mM Tris-HCl, pH 7.5 + 150 mM NaCl + 0.25 mM DTT + 0.1 mM EGTA + 0.1 mM EDTA + 0.1 mM PMSF + 25%		
Purity:	> 95% as determined by densitometry		
Storage & Stability:	Stable for 1 year in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

