

## Recombinant Human CDK8 Inactive GST

<b>Catalog No.</b>	CRC239A	<b>Quantity:</b>	20 µg
	CRC239B		50 µg

**Description:** Recombinant full-length human CDK8 was expressed by baculovirus in Sf9 insect cells using an N-terminal GST tag.

**Concentration:** 0.2 mg/ml

**Protein Accession No:** NM\_005229


**Source:** Sf9 insect cells

**Molecular Weight:** 84 kDa

**Formulation:** Liquid

**Purity:** >80%

**Storage & Stability:** Store product at  $-80^{\circ}\text{C}$ . For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. **Avoid repeated freeze-thaw cycles.**



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.