


SERPINE1

Recombinant Canine PAI-1 Latent Fraction

Catalog No.	CRD004A CRD004B	Quantity:	500 µg 1 mg
Alternate Names:	Plasminogen Activator Inhibitor type 1, SERPINE1, PAI, PAI1, PLANH1		
Description:	Glycosylated wild type dog PAI-1 produced in insect cells. Latent form.		
Concentration:	1 mg/ml		
Gene ID:	403476		
Source:	Insect cells		
Molecular Weight:	46.0 kDa		
Formulation:	Frozen liquid containing 0.05 M Sodium Acetate + 0.1 M NaCl + 1 mM EDTA, pH 5.0		
Purity:	>95% by SDS-PAGE analysis		
Application Note:	Optimal experimental conditions and working dilutions must be determined by the end user.		
Storage & Stability:	Store in working aliquots at -80°C. Stable for 3 years from delivery. Avoid repeated freeze-thaw cycles.		

1. Glycosylated Canine PAI-1 Latent Fraction (3 µg) reduced.
2. Canine PAI-1 (3 µg) + Human µPA (12 µg) reduced.
3. Molecular Weight Markers.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

