


Recombinant Human FABP5 (aa 1-135)

Catalog No.	CRF148A CRF148B	Quantity:	100 µg 500 µg
Alternate Names:	Fatty Acid Binding Protein 5, E-FABP, EFABP, PA-FABP, PAFABP		
Description:	Recombinant Human Fatty Acid Binding Protein 5 containing amino acids 1-135.		
Concentration:	1.0 mg/ml		
Gene ID:	2171		
Protein Accession No:	NP_001435		
Source:	<i>E. coli</i>		
Molecular Weight:	15 kDa		
Formulation:	Liquid in 20 mM Tris, pH 8.0 + 1 mM DTT + 20% glycerol		
Purity:	> 90% as determined by SDS-PAGE analysis		
Amino Acid Sequence:	MATVQQLAGR WRLVDSKGF D EYMKELGVGI ALRKM GAMAK PDCIITCDGK NLTIKTESTL KTTQF SCTLG EKFEETTADG RKTQTVCNFT DGALVQHQEW DGKESTITRK LKDGKLVVEC VMNNVTCTRI YEKVE		
Storage & Stability:	Store in working aliquots at -20°C. Avoid repeated freeze-thaw cycles.		

15% SDS-PAGE

(kDa)



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

