


## Recombinant Human c-JUN (aa 1-79) GST

<b>Catalog No.</b>	CRJ101A CRJ101B	<b>Quantity:</b>	20 µg 50 µg
<b>Description:</b>	Recombinant human c-JUN (1-79) was expressed in E.coli cells using a N-terminal GST tag.		
<b>Concentration:</b>	0.2 mg/ml		
<b>Protein Accession No:</b>	NM_002228		
<b>Source:</b>	<i>E. coli</i>		
<b>Molecular Weight:</b>	35 kDa		
<b>Formulation:</b>	Liquid		
<b>Purity:</b>	>50%		
<b>Storage &amp; Stability:</b>	Store product at $-80^{\circ}\text{C}$ . For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. <b>Avoid repeated freeze-thaw cycles.</b>		



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

