

Recombinant Human PI3K gamma, p110gamma with N-terminal GST-tag Inactive Mutant

Catalog No. CRP059-5 **Quantity:** 5 µg
CRP059 10 µg

Description: The PI3K γ D946A protein is a catalytically inactive mutant of PI3K γ with a D946A mutation in the ATP binding site. This recombinant catalytically inactive protein can be used as a negative control in any kind of PI3K γ kinase activity studies. Recombinant full length PI3K γ DA mutant carries a N terminal GST-tag and was purified by affinity chromatography. The GST-tag facilitates the protein's application in typical GST pull-down assays. Kinase has a MW = 126.3 kDa (without Tag).
Background: Phosphoinositide 3-kinases (PI3Ks) phosphorylate phosphatidylinositols (PIs) at their 3' OH position generating lipid second messengers and thereby regulate numerous biological processes including cell growth, differentiation, survival, proliferation, migration and metabolism. On the basis of structural similarities and substrate specificity, the PI3K family can be subdivided into three classes termed I, II, and III.
All human class I members are heterodimers consisting of a catalytic subunit (MW approx. 110 kDa) and a non-catalytic subunit (MW 50, 55, 85, or 101 kDa) and are known to phosphorylate phosphatidylinositol (PI), phosphatidylinositol-4-mono-phosphate (PIP) and phosphatidylinositol-4,5-bisphosphate (PIP2) in vitro. The class I members can be further subdivided into class IA and IB PI3Ks. Class IA exists in three isoforms (p110 α , p110 β and p110 δ whereas the only class IB member is termed p110 γ . Class IA PI3Ks are activated by adaptor proteins such as Ras or BCAP, or tyrosine-kinase-associated receptors including antigen, co-stimulatory and cytokine receptors (e.g. CD19, CD28, Insulin receptor, EGFR, and PDGFR). p110 γ is activated by G-protein-coupled receptors (GPCRs). Effectors of class I PI3Ks are pleckstrin-homology domain proteins such as Akt/PKB, BTK, TEC, ITK, BAM32, and small GTPases (e.g. Cdc42, Rac, or Ras).
The action of PI3Ks is regulated by the phosphatidyl-inositol-3,4,5-trisphosphate phosphatases SHIP and PTEN.

Source: Sf9 insect cells

Formulation: Liquid. Supplied in 25 mM Tris pH 8.0 + 50 mM NaCl + 0.5 mM MgCl₂ + 50% glycerol.



Cell Sciences, Inc.
480 Neponset Street
Bldg 12A
Canton, MA 02021

Toll Free: 888-769-1246
Phone: 781-828-0610
Fax: 781-828-0542

E-mail: info@cellsciences.com
Web site: www.cellsciences.com

Purity: ≥ 95% by SDS-PAGE.

Storage & Stability: Store at -20°C. **Avoid repeated freeze-thaw cycles.**

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences, Inc.
480 Neponset Street
Bldg 12A
Canton, MA 02021

Toll Free: 888-769-1246
Phone: 781-828-0610
Fax: 781-828-0542

E-mail: info@cellsciences.com
Web site: www.cellsciences.com