

Recombinant Human PAI-1 NBD Mutant

Catalog No.	CRP199B CRP199C	Quantity:	500 µg 1 mg
Alternate Names:	Plasminogen Activator Inhibitor type 1, SERPINE1, PAI, PAI1, PLANH1		
Description:	Recombinant human Plasminogen Activator Inhibitor type 1 NBD labeled mutant contains a Ser to Cys replacement around the vitronectin binding site. Incorporation of an NBD label allows for the reporting of binding to vitronectin.		
Concentration:	2.5 mg/ml		
Gene ID:	5054		
Protein Accession No:	NP_000593		
Source:	<i>E. coli</i>		
Molecular Weight:	43.0 kDa		
Formulation:	Frozen liquid in 0.05 M sodium phosphate + 0.1 M NaCl + 1 mM EDTA, pH 6.6		
Storage & Stability:	Store in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

