

## Recombinant Mouse PAI-1 Latent Fraction

<b>Catalog No.</b>	CRP248B CRP248C	<b>Quantity:</b>	500 µg 1 mg
<b>Alternate Names:</b>	Plasminogen Activator Inhibitor type 1, SERPINE1, PAI, PAI1, PLANH1		
<b>Description:</b>	Recombinant Mouse Plasminogen Activator Inhibitor type 1 Latent Fraction is a useful control for experiments when used in conjunction with the active fraction.		
<b>Concentration:</b>	2.5 mg/ml		
<b>Gene ID:</b>	18787		
<b>Protein Accession No:</b>	NP_032897		
<b>Source:</b>	<i>E. coli</i>		
<b>Molecular Weight:</b>	43.0 kDa		
<b>Formulation:</b>	Liquid in 0.1 M HEPES + 0.1 M NaCl, pH 7.4		
<b>Storage &amp; Stability:</b>	Store in working aliquots at -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

