

## Recombinant Rat PAI-1 Stable Mutant

<b>Catalog No.</b>	CRP252B CRP252C	<b>Quantity:</b>	500 µg 1 mg
<b>Alternate Names:</b>	Plasminogen Activator Inhibitor type 1, SERPINE1, PAI, PAI1, PLANH1		
<b>Description:</b>	Recombinant Rat Plasminogen Activator Inhibitor type 1 Stable Mutant has a single conservative mutation (Isoleucine 91 to Leucine) that yields a rat PAI-1 with an enhanced half-life as compared to the wild type.		
<b>Concentration:</b>	2.5 mg/ml		
<b>Gene ID:</b>	24617		
<b>Protein Accession No:</b>	NP_036752		
<b>Source:</b>	<i>E. coli</i>		
<b>Molecular Weight:</b>	43.0 kDa		
<b>Formulation:</b>	Liquid in 0.05 M sodium acetate + 0.4 M NaCl + 1 mM EDTA, pH 5.0		
<b>Storage &amp; Stability:</b>	Store in working aliquots at -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

