

Recombinant Rat PAI-1 Biotin Latent Fraction

Catalog No.	CRP253B CRP253C	Quantity:	500 µg 1 mg
Alternate Names:	Plasminogen Activator Inhibitor type 1, SERPINE1, PAI, PAI1, PLANH1		
Description:	Recombinant Rat Plasminogen Activator Inhibitor type 1 Biotin Latent Fraction is a useful control for experiments involving the stable mutant form.		
Concentration:	2.5 mg/ml		
Gene ID:	24617		
Protein Accession No:	NP_036752		
Source:	<i>E. coli</i>		
Molecular Weight:	43.0 kDa		
Formulation:	Liquid in 0.1 M HEPES + 0.1 M NaCl, pH 7.4		
Storage & Stability:	Store in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

