



Recombinant Human MEKK2 GST

Catalog No.	CRS164A CRS164B	Quantity:	5 µg 10 µg
Alternate Names:	MAP3K2, MEKK2B		
Description:	Recombinant full-length human MEKK2 was expressed by baculovirus in Sf9 insect cells using an N-terminal GST tag.		
Concentration:	0.1 mg/ml		
Protein Accession No:	NM_006609		
Source:	Sf9 insect cells		
Molecular Weight:	115 kDa		
Formulation:	Liquid in 50 mM Tris-HCl, pH 7.5 + 150 mM NaCl + 0.25 mM DTT + 0.1 mM EGTA + 0.1 mM EDTA + 0.1 mM PMSF + 25% glycerol.		
Purity:	>80% by densitometry		
Specific Activity:	226 nmol /min/mg		
Storage & Stability:	Centrifuge vial prior to opening. Store product at -80°C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.		



Purity: >80% by densitometry. Approx. MW ~115kDa.



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

