

Recombinant Mouse Urokinase N-terminal Fragment

Catalog No.	CRU122B CRU122C	Quantity:	100 µg 1 mg
Alternate Names:	Urokinase Plasminogen Activator, PLAU, ATF, UPA, URK, u-PA		
Description:	N-terminal fragment of recombinant mouse urokinase		
Concentration:	1.0 mg/ml		
Gene ID:	18792		
Protein Accession No:	NP_032899		
Source:	Insect cells		
Molecular Weight:	18.5 kDa		
Formulation:	Liquid in 0.1 M Tris-HCl + 0.15 M NaCl, pH 7.4		
Storage & Stability:	Store in working aliquots at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

